

**POLYTRON® PT 2500 E**



**Stand Dispersing Device (Ecoline)**

**DISPERSING DEVICES FOR BUDGET-CONSCIOUS LABS**

For a variety of applications

« Process volumes from 0.05 to 2,500 ml »



**KINEMATICA**  
DISPERSING AND MIXING TECHNOLOGY

# POLYTRON® - Efficient and Economical

The sleek, sound-insulated casing conceals enough power to safely and quickly handle a variety of applications. Our broad range of dispersing aggregates cover the range from preparing samples to testing product formulas and recipes.

## TECHNICAL INFORMATION

**Processing volume (water)**  
0.05 up to 2,500 ml

**Tip speed**  
Up to max. 26 m/s

**Applicable aggregates**  
ø 3 mm to ø 25 mm

**Speed range**  
500 to 30,000 rpm

**Drive coupling**  
Aggregates with E coupling

**Noise level (drive without aggregate)**  
66 dB(A) at 25,000 rpm  
72 dB(A) at 30,000 rpm

**Motor**  
Universal with carbon brushes

**Drive power**  
500W

**Connection voltage**  
90 - 230 V ± 10%, 50Hz / 60Hz

**Relative humidity (max.)**  
80% during storage  
80% during operation

**Operating temperature**  
0 - 40 °C

**Protection class acc. to DIN**  
IP 20

**Dimensions (L x W x H)**  
Drive: 157 x 76 x 234 mm

**Weight**  
Drive: 1800 g

**EMC standards**  
IEC/EN 61000-6-2 / EN 61000-6-3

**Safety standard**  
IEC/EN 61010-2-51


## THE EXPERT FOR YOUR LABORATORY

For more than 60 years, we have been perfecting the rotor/stator technology, working in close cooperation with users and universities. When using our technology, your test samples will be dispersed into homogeneous, chemical systems. Depending on the type of samples you will get homogeneous emulsions, dispersions or suspensions.

This is the best foundation for all subsequent analyses or formulations.

## ERGONOMICS:

- Easy operation with an on-off switch and speed control
- Intelligent variable speed control for long equipment life
- Speed stability even with changes in viscosity
- Easy adjustment of motor height
- Switch aggregate quickly with Click & Go coupling
- Easy aggregate clean-up with EasyClean Design



## RESULTS ARE CONVINCING

**SMALLEST SIZES OF PARTICLES/DROPLETS**  
because precise analysis require homogeneous samples

**EFFICIENT SAMPLE PREPARATION**  
because fast dispersion saves time

**STANDARD OR CUSTOM-MADE AGGREGATES**  
perfected for and in cooperation with the users

**CLEANING AND DISASSEMBLY**  
EasyClean Design makes it very simple


**QUICK ENGAGEMENT OF THE AGGREGATES** by using Click & Go

**BEST MATERIALS AVAILABLE**  
because without them, we could not call it Swiss Made

**DURABILITY** is a Swiss virtue



# Users / Application ranges



General homogenizer applications

---

Dispersions, suspensions, and emulsions

---

Dispersion of resins and pigments

---

Dispersion of plant, animal and human tissue samples in various liquids and volumes

---

Dissolving solids

---

Preparation of samples for quality assurance/control

---

Breakdown of animal and plant cells

---

Formulation of lotions, creams and foodstuffs

---

Extraction of active pharmaceutical ingredients (API) in pills and tablets

*This list provides a small selection of possible applications.*

# Kinematic dispersion using the rotor/stator system

## Who is the founder of it?

The roots of today's KINEMATICA began in the former chemical/physical research institute of Prof. P. Willems in Lucerne Switzerland, the inventor of the modern rotor/stator dispersion technology.


### THE PRINCIPLE:

The spinning rotor generates a vacuum. The sample is drawn in and discharged to the outside through the stator's slots. Between the rotor/stator (shear gap), the product is subject to high deceleration-tangential and radial acceleration forces. The individual particles/droplets are torn apart and thus reduced in size; they are reduced several micrometers in size through additional cutting and crashing effects.

### THE RESULT:

Microscopic dispersions, emulsions, suspensions, and foams are generated. Droplets, particles, and gas bubbles are reduced to several micrometers or smaller. Small and large substances are dispersed more economically, faster, and better than with any other system.

Your samples are efficiently «**POLYTRONIZED**».



# Aggregate Technology

## Unrivalled!

The high efficiency and the uniqueness of our POLYTRON® dispersing aggregate can only be achieved through superior manufacturing precision and the sophisticated KINEMATICA design method. This is 100% Swiss Made.

### STANDARD AGGREGATES (EC DESIGN WITH TWO GEAR RINGS)

This is all you need for popular applications in the preparation of samples and formulation development. The cutting teeth embedded in the head of the tool optimize the pre-crushing of the sample. An enclosed ring design provides optimum protection for those teeth and prevents them from bending.

### SPECIAL AGGREGATES (W, M, B, G, SYN, X DESIGN)

Specializations in the chemical, pharmaceutical, cosmetic, and food industry, as well as the life sciences, makes new designs inevitable. These aggregates are based on state-of-the-art technology and only available from KINEMATICA. The highly specialized tool is subject to further development that involves close cooperation with our customers and universities. Thus, it is continuously perfected for specific applications.

### MATERIALS

Whether you need a tool with standard or special geometry, all aggregates are manufactured from high-alloy stainless steel 316L, electrolytically polished and made in Switzerland. Our synthetic aggregates (SYN) are manufactured from the highest quality materials available. This guarantees their application at a motor speed up to 30,000 rpm.

### BEARINGS

All dispersing aggregates for the PT 2500 E are equipped with high-quality synthetic bearings. They exhibit excellent mechanical properties and chemical resistance. They are all FDA approved.

### CLEANING AND DISASSEMBLY

EasyClean: Our technology ensures the simplest removal and fast cleaning, either by hand, dishwasher, or autoclave. Disassembling the precisely fitting components is simple and fast.

### LOCKING SYSTEM

All our POLYTRON® aggregates are fitted with a high-precision quick-release system: The Click & Go

# Dispersing aggregates for PT 2500 E

## Your applications, our response

KINEMATICA aggregates solve all problems you are faced with every day in the laboratory. Using our Click & Go quick release makes the changing of aggregates very easy.



### EC «STANDARD DISPERSING AGGREGATE»



#### Design

The universal geometry uses two rows of teeth, is suitable for a large spectrum of laboratory applications, and it is reliable and efficient

#### Selection and applications

- All standard dispersing tasks
- Grinding of animal and human tissue samples in a variety of fluids and volumes
- Sample preparation for the extraction and dissolution of organic materials

- Sample preparation for medical diagnostic procedures as well as DNA analysis
- Extractions of systems/substances from vegetable-derived samples
- Disruption of animal and human cells
- Deagglomeration of united cell structures
- Intensive mixing

Sales No.:	Description	ø/length mm	Volume ml	max. Tip speed m/s
11030004	PT-DA 03/2EC-E050	3 / 50	0.05 - 2	2
11030012	PT-DA 05/2EC-E085	5 / 85	0.1 - 5	5
11030031	PT-DA 07/2EC-E107	7 / 107	0.3 - 10	8
11030062	PT-DA 12/2EC-E123	12 / 123	2 - 250	14
11030042	PT-DA 12/2EC-E157	12 / 157	2 - 250	14
11030069	PT-DA 20/2EC-E192	20 / 192	100 - 2'000	21
11030264	PT-DA 25/2EC-E192	25 / 192	100 - 2'500	26
11030265	PT-DA 25/2FEC-E192*	25 / 192	100 - 2'500	26

\* F Model: For medium to fine particle/droplet sizes

### SYN «SYNTHETICS»



#### Design

Disposable aggregates for the processing of multiple samples per day, where the required sterilization of stainless steel aggregates represents a significant time limitation. Available in PES / Special POM. Sterile packaging in packet sizes of 25 units each (ø7 mm) and 10 units per packet (ø12 mm) each

#### Selection and applications

- Prevention of X contamination
- All standard dispersing applications similar to EC aggregate

Sales No.:	Description	ø/length mm	Volume ml	max. Tip speed m/s
11030030	PT-DA 07/2SYN-E082	7 / 82	0.3 - 10	8
11030035	PT-DA 12/2SYN-E123	12 / 123	2 - 250	14

### X «X-DESIGN»



#### Design

A special geometry designed for dispersing tablets and pills or for preventing suppositories from clogging or agglutinating

#### Selection and applications

- Ultrafast crushing of tablets, coated pills, and suppositories.
- Basis for the analysis of active pharmaceutical ingredients (APIs)
- Used in the substance analysis or for the quality control during the production of tablets

Sales No.:	Description	ø/length mm	Volume ml	max. Tip speed m/s
11030255	PT-DA 20/2XEC-E116	20 / 116	min. 10 *	21

\* When using a vessel ø 22 mm/volume is always dependent on the diameter of the vessel

### M «OUTSIDE CUTTING BLADES»

#### Design

Outside blades pre-crush samples that are larger than the rotors diameter. Subsequently, they can be dispersed inside the generator

#### Selection and applications

- Dispersing of large samples made simple
- All standard dispersing applications similar to EC aggregate



Sales No.:	Description	ø/length mm	Volume ml	max. Tip speed m/s
11030266	PT-DA 12/2MEC-E157	12 / 157	2 - 250	14
11030263	PT-DA 20/2MEC-E192	20 / 192	100 - 2'000	21
11030293	PT-DA 25/2MEC-E192	25 / 192	100 - 2'500	26

### G «GAS TIGHT»

#### Design

Integrated mechanical seals

#### Selection and applications

- Under pressure or vacuum
- Preventing air from being added to emulsions
- Aggressive acid mixtures



Sales No.:	Description	ø/length mm	Volume ml	max. Tip speed m/s
11030273	PT-DA 20/2G-E210	20 / 210	100 - 2'000	21

### W «W-DESIGN»

#### Design

The W design prevents fibrous, stringy or solid samples from clogging the aggregate

#### Selection and applications

- Dispersing fibrous, stringy or solid samples (e.g. pieces of meat)
- All standard dispersing applications similar to EC aggregate



Sales No.:	Description	ø/length mm	Volume ml	max. Tip speed m/s
11030060	PT-DA 12/2WEC-E157	12 / 157	2 - 250	14
11030258	PT-DA 20/2WEC-E192	20 / 192	100 - 2'000	21
11030290	PT-DA 20/2WG-E210*	20 / 210	100 - 2'000	21
11030260	PT-DA 25/2WEC-E192	25 / 192	100 - 2'500	26

\* W-design aggregate in gas-tight model

### B «BIOTRONA® HIGH-TURBULENCE MIXER»

#### Design

High-turbulence mixing head with the least possible shear forces and energy inputs

#### Selection and applications

- Quick removal and suspension of solid particles; also suitable for high viscosity applications



Sales No.:	Description	ø/length mm	Volume ml	max. Tip speed m/s
11030034	PT-DA 07/BEC-E107	7 / 107	10 - 100	8
11030286	PT-DA 12/BEC-E157	12 / 157	50 - 500	14
11030288	PT-DA 20/BEC-E192	20 / 192	100 - 2'000	21
11030289	PT-DA 20/BG-E210*	20 / 210	100 - 2'000	21
11030291	PT-DA 25/BEC-E192	25 / 192	200 - 2'500	26
11030292	PT-DA 25/BG-E210*	25 / 210	200 - 2'500	26

\* BIOTRONA® aggregate in gas-tight model

# Dispersing vessels

Choosing the right dispersing aggregate is crucial for good results. The right vessel also helps contribute to achieving the desired results.

## THE PROBLEM

For normal cylindrical beakers, the flow creates a vortex. The mixing and dispersing efficiency decreases considerably since the product is no longer being optimally processed through the dispersing generator head. This lengthens processing times, requiring more energy input. Depending on the shape of the vortex, this may also lead to additional air being introduced.

## THE SOLUTION

Our specially designed GS jars have a cloverleaf-like shape. This effectively interrupts the flow in the radial direction and generates additional counteracting forces. These high turbulences maximize the mixing and dispersion efficiency, thus minimizing processing time and requiring less energy.

## THE SELECTION

Our GS-glasses are available in chemical- and temperature-resistant borosilicate glass or in stainless steel, with working volumes from a few milliliters to several liters, with or without lid and with or without sealed feedthrough.



GS vessel principle  
(viewed from above)



## THE RIGHT DISPERSING VESSEL

Designation	Sales no.:	Ø / H (mm)	Process volume	Model	Designation	Sales no.:	Ø / H (mm)	Process volume	Model
GS 15	11051010	30 / 100	35 ml	Cloverleaf shaped Borosilicate glass	GS100-1	11052001	105 / 150	1'000 ml	Cylindrical Stainless steel
GS 25	11051011	54 / 150	150 ml		GS100-2	11052002	135 / 170	2'000 ml	
GS 40	11051012	80 / 200	500 ml		GS100-3	11052003	165 / 180	3'000 ml	
GS 50	11051013	90 / 250	1000 ml		GS100-1.C	11052101	105 / 150	1'000 ml	Cylindrical-shaped stainless steel vessel with hook grips and a cover with a drop handle made of stainless steel
GS 60	11051014	100 / 300	1300 ml		GS100-2.C	11052102	135 / 170	2'000 ml	
					GS100-3.C	11052103	165 / 180	3'000 ml	
GS 15 K	11051020	30 / 100	35 ml	Cloverleaf shaped Borosilicate glass	GS100-1.C.B.	11052201	105 / 150	1'000 ml	Cylindrical-shaped stainless steel vessel with hook grips, equipped with a baffles part shelf as an insert as a cover with drop handle
GS 25 K	11051021	54 / 150	150 ml		GS100-2.C.B.	11052202	135 / 170	2'000 ml	
GS 40 K	11051022	80 / 200	500 ml		GS100-3.C.B.	11052203	165 / 180	3'000 ml	
GS 50 K	11051023	90 / 250	1000 ml						
GS 60 K	11051024	100 / 300	1300 ml						
GS 15KL08	11051030	30 / 100	35 ml						
GS 15KL12	11051031	30 / 100	35 ml						
GS 25KL12	11051040	54 / 150	150 ml						
GS 25KL20	11051041	54 / 150	150 ml						
GS 40KL20	11051050	80 / 200	500 ml						
GS 40KL26	11051051	80 / 200	500 ml						
GS 50KL20	11051060	90 / 250	1000 ml						
GS 50KL26	11051061	90 / 250	1000 ml						
GS 60KL20	11051070	100 / 300	1300 ml						
GS 60KL26	11051071	100 / 300	1300 ml						

# Order information and accessories

Depending on your specific requirements, it is possible to expand the PT 2500 E drive with our functional accessories. Our expert personnel will be happy to assist you and answer all the questions you might have.

## ADAPTER

Adapters for use with dispersing aggregate with D coupling

Sales No.:	Designation
11035003	Adapter PT-DA / PT-DA (E/D)

### Description

This adapter connects older D-coupling dispersing aggregates such as PT 2100 with this PT 2500 E motor (with E-coupling).



## DRIVES / COMPLETE SYSTEMS

Sales No.:	Figure	Designation
<b>Motor</b>		
11010055	1	POLYTRON® PT 2500 E, <b>230 V / EU</b> unit cable
11010056	1	POLYTRON® PT 2500 E, <b>230 V / CH</b> unit cable
11010057	1	POLYTRON® PT 2500 E, <b>230 V / UK</b> unit cable
11010058	1	POLYTRON® PT 2500 E, <b>100-120 V /</b> unit cable
<b>Motor is sold with stand and accessories</b>		
11090101	2	PT 2500 E, <b>230 V, EU</b> unit cable, stand ST-P11/600 with vessel holder and safety ring
11090102	2	PT 2500 E, <b>230 V, CH</b> unit cable, stand ST-P11/600 with vessel holder and safety ring
11090103	2	PT 2500 E, <b>230 V, GB</b> unit cable, stand ST-P11/600 with vessel holder and safety ring
11090104	2	PT 2500 E, <b>100-120 V,</b> unit cable, stand ST-P11/600 with vessel holder and safety ring



1



2

## STANDS

Figure	1	2	2	2	2
Sales no.:	11040051	11040033	11040036	11040037	11040034
Designation	ST-P20/600	ST-P 14/590+500, 230 V (EU-plug)	ST-P 14/590+500, 230 V (CH-plug)	ST-P 14/590+500, 230 V (UK-plug)	ST-P 14/590+500, 100-120 V

**Description** Plate stand with two support rods. The motor's height can be adjusted. For more security, the motor cannot be adjusted horizontally.

Electric plate stand with lifting telescopic bar for more ease of use. Incl. up/down button and mat to prevent slipping.

The electronics are designed for ambient temperatures of +5°C to +40°C and were tested per EN 60335-1 and UL 962.



1



2

## STAND ACCESSORIES

Sales no.:	Designation	Description
11045011	Vessel holder with boss head clamp for ST-P20/600	
11045031	Safety ring for stand ST-P20/600	
11036005	Vessel holder with guide for ST-P 14/590+500	
11036001	Extension rod	Short extension rod for PT 2500 E. Required for use with the ST-P 14/590+500 telescopic stand

## YOUR APPLICATION IS OUR FOCUS!

KINEMATICA is a leading manufacturer of dispersing and mixing technology for standard and customized applications in the lab, pilot plant and production areas of pharmaceutical, chemical, food, cosmetic and biotech or life science companies worldwide.

Our POLYTRON® batch and MEGATRON® In-line Homogenizers are suitable for many applications:

- Dispersing of non-soluble liquids or solids into liquids to finest emulsions/multiple emulsions or suspensions
- Induction and dispersion of powders into liquids
- Foaming by gas induction into liquids
- Disintegration of tissue samples for preparation in further analysis
- Dispersing of various samples for quality control

We also deliver POLYMIX® Micro Dry Grinding Mills and a variety of POLYMIX® and MICROTRON® Overhead Stirrers and Mixers.

Whatever your application, we are confident that our team of specialists, with over 60 years of experience, will provide the best solution for you.

## YOUR SATISFACTION IS OUR GOAL!



**KINEMATICA**  
DISPERSING AND MIXING TECHNOLOGY



Headquarters  
**KINEMATICA AG**  
Dispensing and Mixing Technology

Luzernerstrasse 147a  
6014 Luzern  
Switzerland

Tel. +41 41 259 65 65  
Fax +41 41 259 65 75  
E-Mail [info@kinematica.ch](mailto:info@kinematica.ch)

Subsidiaries  
**KINEMATICA, INC.**  
Dispensing and Mixing Technology

155 Keyland Court  
Bohemia, NY 11716  
USA

Tel. +1 631 750 66 53  
Fax +1 631 750 66 57  
E-Mail [info@kinematica-inc.com](mailto:info@kinematica-inc.com)

**KINEMATICA AG**  
Dispensing and Mixing Technology

Office Eastern Europe  
Prosp. Akad. Palladina 44  
UA – Kyiv 03680

Tel./Fax +38 044 422 61 27  
E-Mail [info@kinematica.com.ua](mailto:info@kinematica.com.ua)  
web [www.kinematica.com.ua](http://www.kinematica.com.ua)

**KINEMATICA GMBH**  
Dispensing and Mixing Technology

Münstertäler Str. 12  
79427 Eschbach  
Germany

Tel. +49 7634 504 800 0  
Fax +49 7634 504 800 22  
E-Mail [info@kinematica.de](mailto:info@kinematica.de)

## BRAND NAMES WITH HISTORY



Batch dispersing/homogenizing units



High turbulent mixing



Inline dispersing/homogenizing units



Dispersing/homogenizing Reactors



Stirrers and dry Mills



Mixer/Blender



Viscometers