

Lyovapor L-200 and L-300 Accessories Flyer

Optimize the performance of your Lyovapor L-200 or L-300 with neat accessories.


Beaker Flasks

Simple handling of freeze dryer samples with the beaker flasks. Remove the rubber cover to reach your sample conveniently and maximize the yield. The ideal accessories for all freeze dryers with the manifold set up.

	<p>11066140</p>	<p>Beaker flask 100 ml with rubber cap for direct connection to manifold valves including filter paper (75 µm, 100 pcs.) and O-ring (2 pcs.)</p>
	<p>11066141</p>	<p>Beaker flask 200 ml with rubber cap for direct connection to manifold valves including filter paper (75 µm, 100 pcs.) and O-ring (2 pcs.)</p>
	<p>11066142</p>	<p>Beaker flask 600 ml with rubber cap for direct connection to manifold valves including filter paper (75 µm, 100 pcs.) and O-ring (2 pcs.)</p>
	<p>11066143</p>	<p>Beaker flask 1200 ml with rubber cap for direct connection to manifold valves including filter paper (75 µm, 100 pcs.) and O-ring (2 pcs.)</p>
	<p>11065727</p>	<p>Rubber cap for beaker flask 100 ml, 200 ml</p>
	<p>11065730</p>	<p>Rubber cap for beaker flask 600 ml, 1200 ml</p>

Freeze Drying of Ampules

Freeze drying of ampules with the ampules adapter is absolutely no problem. Just connect your ampules to the connectors and attach them to the manifold freeze dryer.

	11065725	Distributor adapter for 19 ampules
---	-----------------	------------------------------------

Rubber Adapters for Manifolds

To replace current rubber adapters or to change the joint size of the used round bottom flasks.

	11066144	Rubber adapter for manifold valve SJ 29/32 12 pieces
	11066171	Rubber adapter for manifold valve SJ 24/40 12 pieces